

## DATA SHEET

### Sample Pad

**Product** Glass Fibre Sample Pad

#### Basic Information

Composition	Proprietary fiberglass blend; dipped in proprietary chemical blend
Appearance	White, soft paper sheet

#### Technical Properties

Cat No.	Material	Basis Weight (g/m <sup>2</sup> )	Thickness (µm)	Wicking Rate (s/4cm)	Water Absorption (mg/cm <sup>2</sup> )
SP02-03T	100% Glass Fiber	60	0.38	40 ± 15	37.0
SP02-06T	100% Glass Fiber	95	0.60	45 ± 15	46.0
SP02-12T	100% Glass Fiber	155	1.00	40 ± 15	100.0

**Supply Format** Sheet / Strip of desired dimension

**Storage** Temperature: 1° to 35 ° C Material should be kept away from moisture

**Shelf Life** One year from the date of dispatch, stored under above conditions Please retest if using after one year from the date of dispatch.