

---

## DATA SHEET

---

### Absorbent Pad

**Product** Absorbent Pad

#### Basic Information

Composition	100 % cotton linter
Appearance	White transparent

#### Technical Properties

Cat No.	Material	Basis Weight g/m <sup>2</sup>	Thickness (mm)	Wicking Rate s/4 cm	Water Absorption (mg/cm <sup>2</sup> )
1222	100% Cotton Fiber	295	0.83	20	80
1110	100% Cotton Fiber	295	1.16	12	112
1140	100% Cotton Fiber	800	1.40	55	155

**Supply Format** Sheet / Strip of desired dimension

**Storage** Temperature: 1° to 35 ° C  
Material should be kept away from moisture

**Shelf Life** One year from the date of dispatch, stored under above conditions  
Please retest if using after one year from the date of dispatch.